

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-89

Rev No: 03 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Petridish Stand, Acrylic**

**Valid for Catalog No. P10970**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description/Specification

### 1.1 Description: Petridish Stand, Acrylic

1.1.1 Can hold & secure glass & plastic petridishes.

1.1.2 This stand makes transport safe & easier.

1.1.3 Clear racks give the visibility to check culture during incubation.

### 1.2 Design: Rectangular stand

### 1.3 Performance Specifications:

Code	No. of Places	L x W x H (mm)	Compatible	No. Per Case
<b>P10970</b>	28	252 x 105 x 206	90mm Petridish	1

### 1.4 Material:

1.4.1 Made up of Acrylic, confirm US FDA 21CFR

1.4.2 Free from Heavy metals, It is safe use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

### 1.5 Color: Transparent Natural

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-89

Rev No: 03 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Petridish Stand, Acrylic**

## 1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

## 1.7 Basic Features

**1.7.1 Autoclavability:** NO

**1.7.2 Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.7.3 Manufactured in Aseptic Environment**

**1.7.4 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

### Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

### WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ Micro-89

Rev No: 03 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Petridish Stand, Acrylic**

Manufactured & Marketed By :  
**Abdos Labtech Private Limited**  
Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India  
For Complaints/Feedback Contact Customer Care at above address or  
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com  
Lic. No. : MFG/IVD/2024/000115    MADE IN INDIA

IVD UK CA

CE

Electronic IFU

32°C

Scan to visit website



Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

