

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Screw Cap Vial with EPDM O- Ring, PP

Valid for Cat. No.

P10223,P10227,P10231, P10235, P10239, P10223A, P10223A1, P10225R, P10225Y, P10225B, P10225G, P10225O, P10225P, P10260A, P10260A1, P10229R, P10229Y, P10229B, P10229G, P10229O, P10229P, P10270A, P10232A1, P10233R, P10233Y, P10233B, P10233G, P10233O, P10233P, P10272A, P10236A1, P10237R, P10237Y, P10237B, P10237G, P10237O, P10237P, P10274A, P10240A1, P10241R, P10241Y, P10241B, P10241G, P10241O, P10241P, P10222, P10243, P10244, P10245, P10246, P30223, P30227, P30231, P30235, P30239, P30223A, P30223A1, P30225R, P30225Y, P30225B, P30225G, P30225O, P30225P, P30260A, P30260A1, P30229R, P30229Y, P30229B, P30229G, P30229O, P30229P, P30270A, P30232A1, P30233R, P30233Y, P30233B, P30233G, P30233O, P30233P, P30272A, P30236A1, P30237R, P30237Y, P30237B, P30237G, P30237O, P30237P, P30274A, P30240A1, P30241R, P30241Y, P30241B, P30241G, P30241O, P30241P, P30222, P30243, P30244, P30245, P30246, P40223, P40227, P40231, P40235, P40239, P40223A, P40223A1, P40225R, P40225Y, P40225B, P40225G, P40225O, P40225P, P40260A, P40260A1, P40229R, P40229Y, P40229B, P40229G, P40229O, P40229P, P40270A, P40232A1, P40233R, P40233Y, P40233B, P40233G, P40233O, P40233P, P40272A, P40236A1, P40237R, P40237Y, P40237B, P40237G, P40237O, P40237P, P40274A, P40240A1, P40241R, P40241Y, P40241B, P40241G, P40241O, P40241P, P40222, P40243, P40244, P40245, P40246, P50223, P50227, P50231, P50235, P50239, P50223A, P50223A1, P50225R, P50225Y, P50225B, P50225G, P50225O, P50225P, P50260A, P50260A1, P50229R, P50229Y, P50229B, P50229G, P50229O, P50229P, P50270A, P50232A1, P50233R, P50233Y, P50233B, P50233G, P50233O, P50233P, P50272A, P50236A1, P50237R, P50237Y, P50237B, P50237G, P50237O, P50237P, P50274A, P50240A1, P50241R, P50241Y, P50241B, P50241G, P50241O, P50241P, P50222, P50243, P50244, P50245, P50246, P10224, P10228, P10232, P10236, P10240, P10224A, P10224A1, P10226R, P10226Y, P10226Y, P10226B, P10226G, P10226O, P10226P, P10228A, P10228A1, P10230R, P10230Y, P10230B, P10230B, P10230G, P10230O, P10230P, P10232A, P10270A1, P10234R, P10234Y, P10234B, P10234G, P10234O, P10234P, P10236A, P10272A1, P10238R, P10238Y, P10238B, P10238G, P10238O, P10238P, P10240A, P10274A1, P10242R, P10242Y, P10242B, P10242G, P10242O, P10242P, P10221, P10247, P10248, P10249, P10250, P30224, P30228, P30232, P30236, P30240, P30224A, P30224A1, P30226R, P30226Y, P30226B, P30226G, P30226O, P30226P, P30228A, P30228A1, P30230R, P30230Y, P30230B, P30230G, P30230O, P30230P, P30232A, P30270A1, P30234R, P30234Y, P30234B, P30234G,

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30234O, P30234P, P30236A, P30272A1, P30238R, P30238Y, P30238B, P30238G, P30238O, P30238P, P30240A, P30274A1, P30242R, P30242Y, P30242B, P30242G, P30242O, P30242P, P30221, P30247, P30248, P30249, P30250, P40224, P40228, P40232, P40236, P40240, P40224A, P40224A1, P40226R, P40226Y, P40226B, P40226G, P40226O, P40226P, P40228A, P40228A1, P40230R, P40230Y, P40230B, P40230G, P40230O, P40230P, P40232A, P40270A1, P40234R, P40234Y, P40234B, P40234G, P40234O, P40234P, P40236A, P40272A1, P40238R, P40238Y, P40238B, P40238G, P40238O, P40240A, P40274A1, P40242R, P40242Y, P40242B, P40242G, P40242O, P40242P, P40221, P40247, P40248, P40249, P40250, P50224, P50228, P50232, P50236, P50240, P50224A, P50224A1, P50226R, P50226Y, P50226B, P50226G, P50226O, P50226P, P50228A, P50228A1, P50230R, P50230Y, P50230B, P50230G, P50230O, P50230P, P50232A, P50270A1, P50234R, P50234Y, P50234B, P50234G, P50234O, P50234P, P50236A, P50272A1, P50238R, P50238Y, P50238B, P50238G, P50238O, P50238P, P50240A, P50274A1, P50242R, P50242Y, P50242B, P50242G, P50242O, P50242P, P50221, P50247, P50248, P50249, P50250, P10290

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Disposable Screw Cap Vials, DNase, RNase, Human DNA, ATP, Pyrogen Free, PCR Inhibitors & IVD Certified.

1.2 Design:

1.2.1 Vials: Conical / Skirted/Self Stand

1.2.3 Cap: Separate Cap with threads (PP) along with EPDM O- Ring (Black)

1.2.4 without Printing/Printing

1.3 Performance Specifications:

1.3.1 Vial Volume: 0.5ml, 1.5ml & 2.0ml

1.3.2 Brim Volume: 0.8ml, 1.80ml & 2.15ml

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

1.3.3 OD for 0.5, 1.5 & 2ml -10.20mm

1.3.4 Height for 0.5ml – 43.85mm

1.5ml Conical – 43.10mm, 1.5ml Self Stand – 44.0mm

& 2.0ml Conical – 44.30mm, 2.0ml Self stand – 45.10mm

Nonprinted – Non-assembled

Product code	Product Description	Style	Cap & Tube Colour	Printed	Assembled	Sterile	No./Case
P10223	0.5 ml Screw Cap & Tube	self standing	natural	no	no	no	1000
P10224	0.5 ml Screw Cap & Tube	self standing	natural	no	no	yes	1000
P10227	1.5 ml Screw Cap & Tube	self standing	natural	no	no	no	1000
P10228	1.5 ml Screw Cap & Tube	self standing	natural	no	no	yes	1000
P10231	1.5 ml Screw Cap & Tube	conical	natural	no	no	no	1000
P10232	1.5 ml Screw Cap & Tube	conical	natural	no	no	yes	1000
P10235	2.0 ml Screw Cap & Tube	self standing	natural	no	no	no	1000
P10236	2.0 ml Screw Cap & Tube	self standing	natural	no	no	yes	1000
P10239	2.0 ml Screw Cap & Tube	conical	natural	no	no	no	1000
P10240	2.0 ml Screw Cap & Tube	conical	natural	no	no	yes	1000
P10223A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	no	no	no	1000
P10224A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10223A1	0.5 ml Screw Cap & Tube	self standing	cap & tubeTr. Amber	no	no	no	1000
P10224A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	no	yes	1000
P10225R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	no	no	1000
P10225Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	no	no	1000
P10225B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	no	no	1000
P10225G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	no	no	1000
P10225O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	no	no	1000
P10225P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	no	no	1000
P10226R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	no	yes	1000
P10226Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	no	yes	1000
P10226B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	no	yes	1000
P10226G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	no	yes	1000
P10226O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	no	yes	1000
P10226P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	no	yes	1000
P10228A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	no	no	yes	1000
P10228AR	1.5 ml Screw Cap & Tube	self standing	Red cap & tube Amber	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10228AY	1.5 ml Screw Cap & Tube	self standing	Yellow cap & tube Amber	no	no	yes	1000
P10228AB	1.5 ml Screw Cap & Tube	self standing	Blue cap & tube Amber	no	no	yes	1000
P10228AG	1.5 ml Screw Cap & Tube	self standing	Green cap & tube Amber	no	no	yes	1000
P10228AO	1.5 ml Screw Cap & Tube	self standing	Orange cap & tube Amber	no	no	yes	1000
P10228AP	1.5 ml Screw Cap & Tube	self standing	Purple cap & tube Amber	no	no	yes	1000
P10260A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	no	no	no	1000
P10228A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	no	yes	1000
P10260A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	no	no	1000
P10229R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	no	no	1000
P10229Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	no	no	1000
P10229B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	no	no	1000
P10229G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	no	no	1000
P10229O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	no	no	1000
P10229P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	no	no	1000
P10230R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	no	yes	1000
P10230Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10230B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	no	yes	1000
P10230G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	no	yes	1000
P10230O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	no	yes	1000
P10230P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	no	yes	1000
P10232A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	no	no	yes	1000
P10270A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	no	no	no	1000
P10270A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr.Amber	no	no	yes	1000
P10232A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr. Amber	no	no	no	1000
P10233R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	no	no	no	1000
P10233Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	no	no	1000
P10233B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	no	no	no	1000
P10233G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	no	no	no	1000
P10233O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	no	no	no	1000
P10233P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	no	no	no	1000
P10234R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	no	no	yes	1000
P10234Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	no	yes	1000
P10234B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10234G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	no	no	yes	1000
P10234O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	no	no	yes	1000
P10234P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	no	no	yes	1000
P10236A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	no	no	yes	1000
P10272A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	no	no	no	1000
P10272A1	2.0 ml Screw Cap & Tube	self-stand	cap & tubeTr. amber	no	no	yes	1000
P10236A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr.amber	no	no	no	1000
P10237R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	no	no	no	1000
P10237Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	no	no	no	1000
P10237B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	no	no	no	1000
P10237G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	no	no	no	1000
P10237O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	no	no	no	1000
P10237P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	no	no	no	1000
P10238R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	no	no	yes	1000
P10238Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	no	no	yes	1000
P10238B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	no	no	yes	1000
P10238G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10238O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	no	no	yes	1000
P10238P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	no	no	yes	1000
P10240A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	no	no	yes	1000
P10274A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	no	no	no	1000
P10274A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	no	no	yes	1000
P10240A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	no	no	no	1000
P10241R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	no	no	no	1000
P10241Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	no	no	1000
P10241B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	no	no	no	1000
P10241G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	no	no	no	1000
P10241O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	no	no	no	1000
P10241P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	no	no	no	1000
P10242R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	no	no	yes	1000
P10242Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	no	yes	1000
P10242W	2.0 ml Screw Cap & Tube	conical	Cap white & tube natural	no	no	yes	1000
P10242B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	no	no	yes	1000
P10242G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	no	no	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10242O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	no	no	yes	1000
P10242P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	no	no	yes	1000
P10222	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	no	no	no	1000
P10243	1.5ml Screw Cap & Tube LD	Self Stand	Cap & tube natural	no	no	no	1000
P10244	1.5ml Screw Cap & Tube LD	conical	Cap & tube natural	no	no	no	1000
P10245	2.0 ml Screw Cap & Tube LD	Self Stand	Cap & tube natural	no	no	no	1000
P10246	2.0 ml Screw Cap & Tube LD	Conical	Cap & tube natural	no	no	no	1000
P10221	0.5 ml Screw cap & tube LD	Self Stand	Cap & tube natural	no	no	yes	1000
P10247	1.5ml Screw Cap & Tube LD	Self Stand	Cap & tube natural	no	no	yes	1000
P10248	1.5ml Screw Cap & Tube LD	conical	Cap & tube natural	no	no	yes	1000
P10249	2.0 ml Screw Cap & Tube LD	Self Stand	Cap & tube natural	no	no	yes	1000
P10250	2.0 ml Screw Cap & Tube LD	Conical	Cap & tube natural	no	no	yes	1000
P10290	1.5 ml Screw Cap & Tube	Self Stand	Cap & tube natural	no	no	yes	500
P10290R	1.5 ml Screw Cap & Tube	Self Stand	Cap Red & tube natural	no	no	yes	500
P10290Y	1.5 ml Screw Cap & Tube	Self Stand	Cap Yellow & tube natural	no	no	yes	500
P10290G	1.5 ml Screw Cap & Tube	Self Stand	Cap Green & tube natural	no	no	yes	500
P10290B	1.5 ml Screw Cap & Tube	Self Stand	Cap Blue & tube natural	no	no	yes	500

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10290W	1.5 ml Screw Cap & Tube	Self Stand	Cap White & tube natural	no	no	yes	500
P10290O	1.5 ml Screw Cap & Tube	Self Stand	Cap Orange & tube natural	no	no	yes	500
P10290P	1.5 ml Screw Cap & Tube	Self Stand	Cap Purple & tube natural	no	no	yes	500

Printed-Nonassembled

Product code	Product Description	Style	Cap & Tube Colour	Printed	Assembled	No./Case	Sterile
P30223	0.5 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	no
P30224	0.5 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	yes
P30227	1.5 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	no
P30228	1.5 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	yes
P30231	1.5 ml Screw Cap & Tube	conical	natural	Yes	no	1000	no
P30232	1.5 ml Screw Cap & Tube	conical	natural	Yes	no	1000	yes
P30235	2.0 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	no
P30236	2.0 ml Screw Cap & Tube	self standing	natural	Yes	no	1000	yes
P30239	2.0 ml Screw Cap & Tube	conical	natural	Yes	no	1000	no
P30240	2.0 ml Screw Cap & Tube	conical	natural	Yes	no	1000	yes
P30223A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	Yes	no	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30224A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	Yes	no	1000	yes
P30223A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	no	1000	no
P30224A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	no	1000	yes
P30225R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	no	1000	no
P30225Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	no	1000	no
P30225B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	no	1000	no
P30225G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	no	1000	no
P30225O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	no	1000	no
P30225P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	no	1000	no
P30226R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	no	1000	yes
P30226Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	no	1000	yes
P30226B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	no	1000	yes
P30226G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	no	1000	yes
P30226O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	no	1000	yes
P30226P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	no	1000	yes
P30228A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	Yes	no	1000	yes
P30260A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	Yes	no	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30228A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	no	1000	yes
P30260A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	no	1000	no
P30229R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	no	1000	no
P30229Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	no	1000	no
P30229B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	no	1000	no
P30229G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	no	1000	no
P30229O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	no	1000	no
P30229P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	no	1000	no
P30230R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	no	1000	yes
P30230Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	no	1000	yes
P30230B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	no	1000	yes
P30230G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	no	1000	yes
P30230O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	no	1000	yes
P30230P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	no	1000	yes
P30232A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	Yes	no	1000	yes
P30270A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	Yes	no	1000	no
P30270A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr.Amber	Yes	no	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30232A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr. Amber	Yes	no	1000	no
P30233R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	no	1000	no
P30233Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	no	1000	no
P30233B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	no	1000	no
P30233G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	no	1000	no
P30233O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	no	1000	no
P30233P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	no	1000	no
P30234R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	no	1000	yes
P30234Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	no	1000	yes
P30234B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	no	1000	yes
P30234G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	no	1000	yes
P30234O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	no	1000	yes
P30234P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	no	1000	yes
P30236A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	Yes	no	1000	yes
P30272A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	Yes	no	1000	no
P30272A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr. amber	Yes	no	1000	yes
P30236A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr. amber	Yes	no	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30237R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	Yes	no	1000	no
P30237Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	Yes	no	1000	no
P30237B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	Yes	no	1000	no
P30237G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	Yes	no	1000	no
P30237O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	Yes	no	1000	no
P30237P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	Yes	no	1000	no
P30238R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	Yes	no	1000	yes
P30238Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	Yes	no	1000	yes
P30238B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	Yes	no	1000	yes
P30238G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	Yes	no	1000	yes
P30238O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	Yes	no	1000	yes
P30238P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	Yes	no	1000	yes
P30240A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	Yes	no	1000	yes
P30274A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	Yes	no	1000	no
P30274A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	Yes	no	1000	yes
P30240A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	Yes	no	1000	no
P30241R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	no	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30241Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	no	1000	no
P30241B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	no	1000	no
P30241G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	no	1000	no
P30241O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	no	1000	no
P30241P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	no	1000	no
P30242R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	no	1000	yes
P30242Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	no	1000	yes
P30242B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	no	1000	yes
P30242G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	no	1000	yes
P30242O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	no	1000	yes
P30242P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	no	1000	yes
P30222	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	Yes	no	1000	no
P30243	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	no	1000	no
P30244	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	Yes	no	1000	no
P30245	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	no	1000	no
P30246	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	Yes	no	1000	no
P30221	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	Yes	no	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P30247	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	no	1000	yes
P30248	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	Yes	no	1000	yes
P30249	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	no	1000	yes
P30250	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	Yes	no	1000	yes

Nonprinted-Assembled

Product code	Product Description	Style	Cap & Tube Colour	Printed	Assembled	No./Case	Sterile
P40223	0.5 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	no
P40224	0.5 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	yes
P40227	1.5 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	no
P40228	1.5 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	yes
P40231	1.5 ml Screw Cap & Tube	conical	natural	no	Yes	1000	no
P40232	1.5 ml Screw Cap & Tube	conical	natural	no	Yes	1000	yes
P40235	2.0 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	no
P40236	2.0 ml Screw Cap & Tube	self standing	natural	no	Yes	1000	yes
P40239	2.0 ml Screw Cap & Tube	conical	natural	no	Yes	1000	no
P40240	2.0 ml Screw Cap & Tube	conical	natural	no	Yes	1000	yes
P40223A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	no	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40224A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	no	Yes	1000	yes
P40223A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	Yes	1000	no
P40224A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	Yes	1000	yes
P40225R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	Yes	1000	no
P40225Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	Yes	1000	no
P40225B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	Yes	1000	no
P40225G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	Yes	1000	no
P40225O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	Yes	1000	no
P40225P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	Yes	1000	no
P40226R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	Yes	1000	yes
P40226Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	Yes	1000	yes
P40226B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	Yes	1000	yes
P40226G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	Yes	1000	yes
P40226O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	Yes	1000	yes
P40226P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	Yes	1000	yes
P40228A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	no	Yes	1000	yes
P40260A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	no	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40228A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	Yes	1000	yes
P40260A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	no	Yes	1000	no
P40229R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	Yes	1000	no
P40229Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	Yes	1000	no
P40229B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	Yes	1000	no
P40229G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	Yes	1000	no
P40229O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	Yes	1000	no
P40229P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	Yes	1000	no
P40230R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	no	Yes	1000	yes
P40230Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	no	Yes	1000	yes
P40230B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	no	Yes	1000	yes
P40230G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	no	Yes	1000	yes
P40230O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	no	Yes	1000	yes
P40230P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	no	Yes	1000	yes
P40232A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	no	Yes	1000	yes
P40270A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	no	Yes	1000	no
P40270A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr.Amber	no	Yes	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40232A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr. Amber	no	Yes	1000	no
P40233R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	no	Yes	1000	no
P40233Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	Yes	1000	no
P40233B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	no	Yes	1000	no
P40233G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	no	Yes	1000	no
P40233O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	no	Yes	1000	no
P40233P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	no	Yes	1000	no
P40234R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	no	Yes	1000	yes
P40234Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	Yes	1000	yes
P40234B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	no	Yes	1000	yes
P40234G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	no	Yes	1000	yes
P40234O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	no	Yes	1000	yes
P40234P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	no	Yes	1000	yes
P40236A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	no	Yes	1000	yes
P40272A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	no	Yes	1000	no
P40272A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr. amber	no	Yes	1000	yes
P40236A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr.amber	no	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40237R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	no	Yes	1000	no
P40237Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	no	Yes	1000	no
P40237B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	no	Yes	1000	no
P40237G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	no	Yes	1000	no
P40237O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	no	Yes	1000	no
P40237P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	no	Yes	1000	no
P40238R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	no	Yes	1000	yes
P40238Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	no	Yes	1000	yes
P40238B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	no	Yes	1000	yes
P40238G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	no	Yes	1000	yes
P40238O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	no	Yes	1000	yes
P40238P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	no	Yes	1000	yes
P40240A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	no	Yes	1000	yes
P40274A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	no	Yes	1000	no
P40274A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	no	Yes	1000	yes
P40240A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	no	Yes	1000	no
P40241R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	no	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40241Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	Yes	1000	no
P40241B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	no	Yes	1000	no
P40241G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	no	Yes	1000	no
P40241O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	no	Yes	1000	no
P40241P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	no	Yes	1000	no
P40242R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	no	Yes	1000	yes
P40242Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	no	Yes	1000	yes
P40242B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	no	Yes	1000	yes
P40242G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	no	Yes	1000	yes
P40242O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	no	Yes	1000	yes
P40242P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	no	Yes	1000	yes
P40222	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	no	Yes	1000	no
P40243	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	no	Yes	1000	no
P40244	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	no	Yes	1000	no
P40245	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	no	Yes	1000	no
P40246	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	no	Yes	1000	no
P40221	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	no	Yes	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P40247	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	no	Yes	1000	yes
P40248	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	no	Yes	1000	yes
P40249	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	no	Yes	1000	yes
P40250	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	no	Yes	1000	yes

Printed-Assembled

Product code	Product Description	Style	Cap & Tube Colour	Printed	Assembled	No./Case	Sterile
P50223	0.5 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	no
P50224	0.5 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	yes
P50227	1.5 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	no
P50228	1.5 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	yes
P50231	1.5 ml Screw Cap & Tube	conical	natural	Yes	Yes	1000	no
P50232	1.5 ml Screw Cap & Tube	conical	natural	Yes	Yes	1000	yes
P50235	2.0 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	no
P50236	2.0 ml Screw Cap & Tube	self standing	natural	Yes	Yes	1000	yes
P50239	2.0 ml Screw Cap & Tube	conical	natural	Yes	Yes	1000	no
P50240	2.0 ml Screw Cap & Tube	conical	natural	Yes	Yes	1000	yes
P50223A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	Yes	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50224A	0.5 ml Screw Cap & Tube	self standing	cap & tube amber	Yes	Yes	1000	yes
P50223A1	0.5 ml Screw Cap & Tube	self standing	cap & tubeTr. Amber	Yes	Yes	1000	no
P50224A1	0.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	Yes	1000	yes
P50225R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	Yes	1000	no
P50225Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	Yes	1000	no
P50225B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	Yes	1000	no
P50225G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	Yes	1000	no
P50225O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	Yes	1000	no
P50225P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	Yes	1000	no
P50226R	0.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	Yes	1000	yes
P50226Y	0.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	Yes	1000	yes
P50226B	0.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	Yes	1000	yes
P50226G	0.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	Yes	1000	yes
P50226O	0.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	Yes	1000	yes
P50226P	0.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	Yes	1000	yes
P50228A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	Yes	Yes	1000	yes
P50260A	1.5 ml Screw Cap & Tube	self standing	cap & tube Amber	Yes	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50228A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	Yes	1000	yes
P50260A1	1.5 ml Screw Cap & Tube	self standing	cap & tube Tr. Amber	Yes	Yes	1000	no
P50229R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	Yes	1000	no
P50229Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	Yes	1000	no
P50229B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	Yes	1000	no
P50229G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	Yes	1000	no
P50229O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	Yes	1000	no
P50229P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	Yes	1000	no
P50230R	1.5 ml Screw Cap & Tube	self standing	cap red & tube natural	Yes	Yes	1000	yes
P50230Y	1.5 ml Screw Cap & Tube	self standing	cap yellow & tube natural	Yes	Yes	1000	yes
P50230B	1.5 ml Screw Cap & Tube	self standing	cap blue & tube natural	Yes	Yes	1000	yes
P50230G	1.5 ml Screw Cap & Tube	self standing	cap green & tube natural	Yes	Yes	1000	yes
P50230O	1.5 ml Screw Cap & Tube	self standing	cap orange & tube natural	Yes	Yes	1000	yes
P50230P	1.5 ml Screw Cap & Tube	self standing	cap purple & tube natural	Yes	Yes	1000	yes
P50232A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	Yes	Yes	1000	yes
P50270A	1.5 ml Screw Cap & Tube	conical	cap & tube Amber	Yes	Yes	1000	no
P50270A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr.Amber	Yes	Yes	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50232A1	1.5 ml Screw Cap & Tube	conical	cap & tube Tr. Amber	Yes	Yes	1000	no
P50233R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	Yes	1000	no
P50233Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	Yes	1000	no
P50233B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	Yes	1000	no
P50233G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	Yes	1000	no
P50233O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	Yes	1000	no
P50233P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	Yes	1000	no
P50234R	1.5 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	Yes	1000	yes
P50234Y	1.5 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	Yes	1000	yes
P50234B	1.5 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	Yes	1000	yes
P50234G	1.5 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	Yes	1000	yes
P50234O	1.5 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	Yes	1000	yes
P50234P	1.5 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	Yes	1000	yes
P50236A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	Yes	Yes	1000	yes
P50272A	2.0 ml Screw Cap & Tube	self-stand	cap & tube amber	Yes	Yes	1000	no
P50272A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr. amber	Yes	Yes	1000	yes
P50236A1	2.0 ml Screw Cap & Tube	self-stand	cap & tube Tr.amber	Yes	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50237R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	Yes	Yes	1000	no
P50237Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	Yes	Yes	1000	no
P50237B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	Yes	Yes	1000	no
P50237G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	Yes	Yes	1000	no
P50237O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	Yes	Yes	1000	no
P50237P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	Yes	Yes	1000	no
P50238R	2.0 ml Screw Cap & Tube	self-stand	cap red & tube natural	Yes	Yes	1000	yes
P50238Y	2.0 ml Screw Cap & Tube	self-stand	cap yellow & tube natural	Yes	Yes	1000	yes
P50238B	2.0 ml Screw Cap & Tube	self-stand	cap blue & tube natural	Yes	Yes	1000	yes
P50238G	2.0 ml Screw Cap & Tube	self-stand	cap green & tube natural	Yes	Yes	1000	yes
P50238O	2.0 ml Screw Cap & Tube	self-stand	cap orange & tube natural	Yes	Yes	1000	yes
P50238P	2.0 ml Screw Cap & Tube	self-stand	cap purple & tube natural	Yes	Yes	1000	yes
P50240A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	Yes	Yes	1000	yes
P50274A	2.0 ml Screw Cap & Tube	conical	cap & tube amber	Yes	Yes	1000	no
P50274A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	Yes	Yes	1000	yes
P50240A1	2.0 ml Screw Cap & Tube	conical	cap & tube Tr. amber	Yes	Yes	1000	no
P50241R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	Yes	1000	no

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50241Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	Yes	1000	no
P50241B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	Yes	1000	no
P50241G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	Yes	1000	no
P50241O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	Yes	1000	no
P50241P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	Yes	1000	no
P50242R	2.0 ml Screw Cap & Tube	conical	cap red & tube natural	Yes	Yes	1000	yes
P50242Y	2.0 ml Screw Cap & Tube	conical	cap yellow & tube natural	Yes	Yes	1000	yes
P50242B	2.0 ml Screw Cap & Tube	conical	cap blue & tube natural	Yes	Yes	1000	yes
P50242G	2.0 ml Screw Cap & Tube	conical	cap green & tube natural	Yes	Yes	1000	yes
P50242O	2.0 ml Screw Cap & Tube	conical	cap orange & tube natural	Yes	Yes	1000	yes
P50242P	2.0 ml Screw Cap & Tube	conical	cap purple & tube natural	Yes	Yes	1000	yes
P50222	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	no
P50243	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	no
P50244	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	Yes	Yes	1000	no
P50245	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	no
P50246	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	Yes	Yes	1000	no
P50221	0.5 ml Screw cap & tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	yes

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P50247	1.5ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	yes
P50248	1.5ml Screw Cap & Tube LD	conical	cap & tube natural	Yes	Yes	1000	yes
P50249	2.0 ml Screw Cap & Tube LD	Self Stand	cap & tube natural	Yes	Yes	1000	yes
P50250	2.0 ml Screw Cap & Tube LD	Conical	cap & tube natural	Yes	Yes	1000	yes

Non Printed Tube with Looped Cap

Product code	Capacity(ml)	Outer Diameter(mm)	Height(mm)	Style	Sterile	Pcs/case
P10252	0.5 ml	10.2	43.85	self standing	no	1000
P10253	0.5 ml	10.2	43.85	self standing	yes	1000
P10254	1.5 ml	10.2	44	self standing	no	1000
P10255	1.5 ml	10.2	44	self standing	yes	1000
P10256	1.5 ml	10.2	43.1	conical	no	1000
P10257	1.5 ml	10.2	43.1	conical	yes	1000
P10258	2.0 ml	10.2	45.1	self standing	no	1000
P10259	2.0 ml	10.2	45.1	self standing	yes	1000
P10260	2.0 ml	10.2	44.3	conical	no	1000
P10261	2.0 ml	10.2	44.3	conical	yes	1000
P10252A	0.5 ml amber	10.2	43.85	self standing	no	1000
P10253A	0.5 ml amber	10.2	43.85	self standing	yes	1000
P10252A1	0.5 ml transparent amber	10.2	43.85	self standing	no	1000
P10253A1	0.5 ml transparent amber	10.2	43.85	self standing	yes	1000
P10262	0.5 ml colour	10.2	43.85	self standing	no	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10263	0.5 ml coloured	10.2	43.85	self standing	yes	1000
P10255A	1.5 ml amber	10.2	43.95	self standing	yes	1000
P10254A	1.5 ml amber	10.2	44	self standing	no	1000
P10255A1	1.5 ml transparent amber	10.2	44	self standing	yes	1000
P10254A1	1.5 ml transparent amber	10.2	44	self standing	no	1000
P10264	1.5 ml colour	10.2	44	self standing	no	1000
P10265	1.5 ml coloured	10.2	44	self standing	yes	1000
P10257A	1.5 ml amber	10.2	43.1	conical	yes	1000
P10256A	1.5 ml amber	10.2	43.1	conical	no	1000
P10257A1	1.5 ml transparent amber	10.2	43.1	conical	yes	1000
P10256A1	1.5 ml transparent amber	10.2	43.1	conical	no	1000
P10266	1.5 ml coloured	10.2	43.1	conical	no	1000
P10267	1.5 ml coloured	10.2	43.1	conical	yes	1000
P10259A	2.0 ml amber	10.2	45.1	self-stand	yes	1000
P10258A	2.0 ml amber	10.2	45.1	self-stand	no	1000
P10259A1	2.0 ml transparent amber	10.2	45.1	self-stand	yes	1000
P10258A1	2.0 ml transparent amber	10.2	45.1	self-stand	no	1000
P10268	2.0 ml coloured	10.2	45.1	self-stand	no	1000
P10269	2.0 ml coloured	10.2	45.1	self-stand	yes	1000
P10261A	2.0 ml amber	10.2	44.3	conical	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10268A	2.0 ml amber	10.2	44.3	conical	no	1000
P10261A1	2.0 ml transparent amber	10.2	44.3	conical	yes	1000
P10268A1	2.0 ml transparent amber	10.2	44.3	conical	no	1000
P10270	2.0 ml coloured	10.2	44.3	conical	no	1000
P10271	2.0 ml coloured	10.2	44.3	conical	yes	1000
P10272	0.5 ml LD	10.2	43.85	Self Stand	no	1000
P10273	1.5ml LD	10.2	44	Self Stand	no	1000
P10274	1.5ml LD	10.2	43.1	conical	no	1000
P10275	2.0 ml LD	10.2	45.1	Self Stand	no	1000
P10276	2.0 ml LD	10.2	44.3	Conical	no	1000
P10277	0.5 ml LD	10.2	43.85	Self Stand	yes	1000
P10278	1.5ml LD	10.2	44	Self Stand	yes	1000
P10279	1.5ml LD	10.2	43.1	conical	yes	1000
P10280	2.0 ml LD	10.2	45.1	Self Stand	yes	1000
P10281	2.0 ml LD	10.2	44.3	Conical	yes	1000

Printed Tube with Looped Cap

Product code	Capacity(ml)	Outer Diameter(mm)	Height (mm)	Style	Sterile	Pcs/case
P10340	0.5 ml	10.2	43.85	self standing	no	1000
P10341	0.5 ml	10.2	43.85	self standing	yes	1000
P10342	1.5 ml	10.2	44	self standing	no	1000
P10343	1.5 ml	10.2	44	self standing	yes	1000
P10344	1.5 ml	10.2	43.1	conical	no	1000
P10345	1.5 ml	10.2	43.1	conical	yes	1000
P10346	2.0 ml	10.2	45.1	self standing	no	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10347	2.0 ml	10.2	45.1	self standing	yes	1000
P10348	2.0 ml	10.2	44.3	conical	no	1000
P10349	2.0 ml	10.2	44.3	conical	yes	1000
P10340A	0.5 ml amber	10.2	43.85	self standing	no	1000
P10341A	0.5 ml amber	10.2	43.85	self standing	yes	1000
P10340A1	0.5 ml transparent amber	10.2	43.85	self standing	no	1000
P10341A1	0.5 ml transparent amber	10.2	43.85	self standing	yes	1000
P10350	0.5 ml coloured	10.2	43.85	self standing	no	1000
P10351	0.5 ml coloured	10.2	43.85	self standing	yes	1000
P10343A	1.5 ml amber	10.2	43.95	self standing	yes	1000
P10342A	1.5 ml amber	10.2	44	self standing	no	1000
P10343A1	1.5 ml transparent amber	10.2	44	self standing	yes	1000
P10342A1	1.5 ml transparent amber	10.2	44	self standing	no	1000
P10352	1.5 ml coloured	10.2	44	self standing	no	1000
P10353	1.5 ml coloured	10.2	44	self standing	yes	1000
P10345A	1.5 ml amber	10.2	43.1	conical	yes	1000
P10344A	1.5 ml amber	10.2	43.1	conical	no	1000
P10345A1	1.5 ml transparent amber	10.2	43.1	conical	yes	1000
P10344A1	1.5 ml transparent amber	10.2	43.1	conical	no	1000
P10354	1.5 ml coloured	10.2	43.1	conical	no	1000
P10355	1.5 ml coloured	10.2	43.1	conical	yes	1000

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

P10347A	2.0 ml amber	10.2	45.1	self-stand	yes	1000
P10346A	2.0 ml amber	10.2	45.1	self-standing	no	1000
P10347A1	2.0 ml transparent amber	10.2	45.1	self-stand	yes	1000
P10346A1	2.0 ml transparent amber	10.2	45.1	self-standing	no	1000
P10356	2.0 ml coloured	10.2	45.1	self-standing	no	1000
P10357	2.0 ml coloured	10.2	45.1	self-stand	yes	1000
P10349A	2.0 ml amber	10.2	44.3	conical	yes	1000
P10348A	2.0 ml amber	10.2	44.3	conical	no	1000
P10349A1	2.0 ml transparent amber	10.2	44.3	conical	yes	1000
P10348A1	2.0 ml transparent amber	10.2	44.3	conical	no	1000
P10358	2.0 ml coloured	10.2	44.3	conical	no	1000
P10359	2.0 ml coloured	10.2	44.3	conical	yes	1000
P10360	0.5 ml LD	10.2	43.85	Self Standing	no	1000
P10361	1.5ml LD	10.2	44	Self Standing	no	1000
P10362	1.5ml LD	10.2	43.1	conical	no	1000
P10363	2.0 ml LD	10.2	45.1	Self Standing	no	1000
P10364	2.0 ml LD	10.2	44.3	Conical	no	1000
P10365	0.5 ml LD	10.2	43.85	Self Stand	yes	1000
P10366	1.5ml LD	10.2	44	Self Stand	yes	1000
P10367	1.5ml LD	10.2	43.1	conical	yes	1000
P10368	2.0 ml LD	10.2	45.1	Self Stand	yes	1000
P10369	2.0 ml LD	10.2	44.3	Conical	yes	1000

1.3.5 Graduation and Writing Surface:

1.3.5.1 Large Printed Graduation

1.3.5.2 White Patch (Frosted) for easy writing on the tube & Chloroform proof

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

1.4 Material:

1.4.1 Polypropylene (PP), Hydrophobic, confirm US FDA 21CFR

1.4.2 The product complies with USP class VI. Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.4.4 This Material is Gamma-Ray sterilizable & Phthalates Free.

1.4.5 EPDM ring meets the FDA criteria Regulation 21 CFR.

1.5 Color:

Natural/ Red/Yellow/Blue/Green/Orange/Purple/Amber/White/Tr. Amber/Last Drop

1.6 Sterilization:

Sterile catalog number from this lot have been pre-sterilizes and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation, SAL 10⁻⁶

1.7 Quality Control:

Products are manufactured, inspected, and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.8 Other information:

Recommended for single use only

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

Pack Size:

Non-Assembled

Vial 100 Pieces x 10pouch + Cap with EPDM O Ring - 100 pieces x 10pouch =1000 per case (20pouch)

Assembled

Vial+Cap 100 Pieces x 10pouch with EPDM O Ring =1000 per case (10pouch)

2. Basic Features

2.1 Temperature range: Tube: -196 °C to + 121 °C

* Caution: Do not use Abdos Cryogenic Vial in the liquid phase of liquid nitrogen. Improper use may cause liquified nitrogen to be trapped inside the vial and lead to pressure build up resulting in possible explosion or biohazard release. Can be used in vapour phase of liquid nitrogen.

2.2 Autoclavability: Autoclavable at 15 psi and 121 °C for 20 minutes (Caps & Vials)

2.3 Leak Test at 635 mm Hg with Cap fitted: No leakage

2.4. These Vials have EPDM 'O' ring to assure efficient sealing and avoid evaporation.

2.5 Centrifugation: 20,000 RCF (* Rated with water)

2.6 Manufactured in Clean room ISO Class 8 (100,000)

2.7 Chemical resistance: see website (www.abdoslifesciences.com)

2.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Autoclaving Guideline for Plastics

TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP


All plastic consumables are Autoclavable at 121⁰C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness, pouch punctured etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimen.
Ensure all sample tubes are evenly filled
For each tube inserted in the rotor, add a tube of equal weight directly opposite it.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Wear a face shield and / or safety goggles if you have to work anywhere near a centrifuge that's in use.
Do not bump, jar, or move the centrifuge while the rotor is spinning. Make sure you don't have the cord dangling from a table edge where someone could catch their foot in it and pull down the centrifuge
Consult local fire codes prior to storage of flammable liquids.



TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Screw Cap Vial with EPDM O- Ring, PP

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110



TECHNICAL DATA SHEET

Doc No – Ab/TDS/CT-03
Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Screw Cap Vial with EPDM O- Ring, PP

