

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-237
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEX – Classic Rotating Mixer

Valid for Catalog No. E11305

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description:

1.1.1 Wavex tube rotators provide gentle but effective mixing of biological samples, held in 1.5 ml to 50 ml tubes.

1.1.2 Used in a variety of applications including immune precipitations, prevention of blood coagulation and latex diagnostics.

1.1.3 Wavex-Classic Tube Rotator has Disk Type Design.

1.2 Design:

Rotisserie Type

1.3 performance Specifications:

1.3.1 Voltage (AC)	: 220~240 Volt
1.3.2 Frequency (Hz)	: 50-60
1.3.3 Power (W)	: 20
1.3.4 Run type	: Continuous operation
1.3.5 Motor type	: DC Motor
1.3.6 Speed range (RPM)	: 0 ~ 80
1.3.7 Overall dimension	: 150 x 530 x 190 mm
1.3.8 Permissible ambient Temp.	: 5- 40°C
1.3.9 Permissible Relative humidity:	80 %

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-237
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEX – Classic Rotating Mixer

1.4.0 Weight	: 8kg
1.4.1 Speed & Time Display: Scale	

1.4 Material: Body – ABS

1.5 Color: Dark Blue

1.6. Spare Accessories:

E11301 Plastic clips for 1.5ml x 12 MCT – 1pcs/case
E11302 Plastic clips for 15ml x 12 MCT – 1pcs/case
E11303 Plastic clips for 50ml x 12 MCT – 1pcs/case

1.6.2 CE certified

1.6.3 **Chemical resistance:** see website (www.abdoslifesciences.com)

1.7 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Storage temp. Max 30°C
Before using the product check crackness, etc.
Check all the information printed on pouches or boxes.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-237
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
WAVEX – Classic Rotating Mixer

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure
Consult local fire codes prior to storage of flammable liquids.



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

