

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-244  
Rev No: 06 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**WAVEMAX – Orbital Shaker**

**Valid for Catalog No. E11290**

**SRN (Single Registration Number): IN-MF-000022272**

## 1. Description Specifications:

### 1.1 Description: WAVEMAX- Orbital Shaker

1.1.1 Abdos powerful orbital shaker is capable of mixing liquids upto a load capacity of 6 kg in bottles, flasks, test tubes & dishes.

1.1.1 Abdos orbital shaker has digital speed control of upto 300 rpm.

1.1.2 These shakers have high capacity of upto 24 x 125 ml, 14 x 250 ml, 9 x 500 ml, 5 x 1l and 2 x 2l flasks and bottles.

1.1.3 These shakers have horizontal circular 19 mm orbit for aeration and mixing.

1.1.4 Universal magnetic clamp can be instantly exchanged.

### 1.3 Performance Specifications:

Specifications	E11290
Voltage	220-240V
Frequency	50 Hz
Motor Type	DC Motor
Shaking Motion	Orbital
Orbit	19 mm
Max Load Capacity (with Platform)	6.4 kg
Speed range	30-300 rpm, variable
Operation Type	LCD
Dimension	280x410x75 mm

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-244  
Rev No: 06 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**WAVEMAX – Orbital Shaker**

Weight	11.25 kg
Permissible ambient Temp.	+4 to +65°C
Permissible ambient Humidity	80%
Speed Increment	1 rpm

**1.4 Material:** Body – ABS

**1.5 Color:** Grey

**1.6 Classic Orbital Shaker**

**E11290** – Wavemax - Orbital Shaker with flat mat - 13x12 inch

**Spare Accessories:**

<b>E11291</b> – Universal magnetic platform for flasks & tube racks - 14x12 inch

**1.7.** CE certified

**1.8 Chemical resistance:** see website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

**1.9 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

**Warranty:** 12 months from the date of Installation or 15 months from the date of Invoicing.

## **Instruction for Use**

Storage temp. Max 30°C

Before using the product check crackness, etc.

Check all the information printed on pouches or boxes.

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/Eqpt-244  
Rev No: 06 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**WAVEMAX – Orbital Shaker**

## WARNING

**DO NOT** store strong oxidizing agents.

**DO NOT** place any plastic labware in a flame.

**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

