

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-162  
Rev No: 08 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Scoop, HDPE**

Valid for Catalog No. P11901, P11902, P11903, P11904, P11905

**SRN (Single Registration Number):** IN-MF-000022272

## 1. Description/Specification

### 1.1 Description: Scoop, HDPE

1.1.1 These scoops used for withdrawal and movement of solids, granulates, gels and other sample.

1.1.2 Round bottom for better use.

1.1.3 It has comfortable, stable and robust handle with precise filling edge.

1.1.4 Natural in colour & superior chemical resistance.

1.1.5 Suitable for use with food.

### 1.2 Design: Round bottom Scoop

### 1.3 Performance Specifications:

Code	Sterile	Capacity (ml)	Total Length (mm)	No. Per Case
P11901	NO	90	180	12
P11902	NO	250	240	12
P11903	NO	600	300	6
P11904	NO	1100	360	6
P11905	NO	1500	400	3
P11901S	YES	90	180	1
P11902S	YES	250	240	1

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-162  
Rev No: 08 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Scoop, HDPE**

<b>P11903S</b>	YES	600	300	1
<b>P11904S</b>	YES	1100	360	1
<b>P11905S</b>	YES	1500	400	1

## 1.4 Material:

1.4.1 High Density poly Ethylene (HDPE), Hydrophobic, confirm US FDA 21CFR

1.4.2 Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

**1.5 Color:** Translucent (Natural)

**1.6 Quality control:** Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release amendments.

## 1.7 Other information

1.7.1 Impact resistant and tough.

1.7.2 Superior for its hardness, gloss, toughness, and electrical insulation properties.

1.7.3 **Chemical resistance:** See website ([www.abdoslifesciences.com](http://www.abdoslifesciences.com))

1.7.4 **Manufactured** in Aseptic Environment

# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-162  
Rev No: 08 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Scoop, HDPE**

**1.8 Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.


**1.9 Sterilization:** Sterile catalog number from this lot have been pre-sterilizes and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products- Requirements for validation and routine control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation, SAL 10<sup>-6</sup>

## Instruction for Use

Product can be used for food contact also.  
Storage temp. Max 30°C  
Before using the product check crackness etc.  
Check all the information printed on pouches or boxes.  
Always wear gloves while handling specimens.

## WARNING

**DO NOT** store strong oxidizing agents.  
**DO NOT** place any plastic labware in a flame.  
**DO NOT** mix any chemicals that may result in a thermic reaction, which can cause product failure  
Consult local fire codes prior to storage of flammable liquids.



# TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-162  
Rev No: 08 Date- 05/2024-25

**Abdos Labtech Pvt. Ltd.**  
**Scoop, HDPE**

---

Quality Manager  
Abdos Labtech Private Limited  
Employee Code 30110

