

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-191

Rev No: 05 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet

Valid for Catalog No.

P80183, P80184, P80185, P80186

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet

1.2 Design :Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet

1.3 Performance Specification:

1.3.1 Best suited for test tubes.

1.3.2 For the stirring of liquids on magnetic stirrers & hot plates stirrers.

1.3.3 Strong alloy magnets sheathed in pure white inert PTFE.

1.4 Material:

1.4.1 Moulded in from FDA approved PTFE & ALNICO v Magnet.

1.5 Color: White

1.6 Quality Control: Testing of various characteristics in accordance with the valid Specification.

1.7 Basic Features:

1.7.1 Made up of Highly Pure Virgin PTFE Material.

1.7.2 Resistance to virtual all chemical attacks.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-191

Rev No: 05 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.

Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet

Code	Type	OD x L (mm)	No. per Case	Basic Features
P80183	Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet	8x11	2	- For the stirring of liquids on magnetic stirrers & hot plates stirrers -Resistance with chemicals & with stand upto 260°C
P80184		10x16	2	
P80185		13x10	2	
P80186		13x20	2	

1.8 Basic Features:

1.8.1 Chemical resistance: see website (www.abdoslifesciences.com)

1.8.2 Manufactured in Aseptic Environment

1.8.3 Autoclavability: Yes

1.8.4 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ GLP-191

Rev No: 05 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.

Wings Magnetic Stirrer Bar, PTFE/ALNICO V Magnet

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

