

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-37
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Fixed Premium Pipettes, PC (Non-Autoclavable)

**Valid for Catalog No. E10001, E10002, E10003, E10004, E10005, E10006, E10007,
E10008, E10009, E10010**

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Fixed Premium Pipettes, PC

- 1.1.1 Abdos fixed premium pipettes are very accurate and gives desired results.
- 1.1.2 Abdos fixed premium pipettes are lightweight, ergonomic, with low force design.
- 1.1.3 Volume settings can be clearly read through digital display.
- 1.1.4 Abdos fixed premium pipettes are easy to calibrate & maintain with the tool supplied.
- 1.1.5 Abdos fixed premium pipettes are designed to avoid repetitive strain injuries.
- 1.1.6 Abdos fixed premium pipettes are supplied with operator manual that provides complete information on calibration & maintenance of the pipettes.
- 1.1.7 Abdos fixed premium pipettes are suitable for transfer liquid volumes.
- 1.1.8 Matching universal style pipette tips, calibrated in accordance with ISO 8655. The quality control involves gravimetric testing of each pipette with distilled water at 22°C.

1.2 Design: Light weight & simple to use operation.

1.3 Performance Specifications: Push Button, Cap, Push Button Rod, Tip ejector, Grip cover, Handle, Display cover, Tip ejector collar, Tip cone, Volume setting, Sealing & ejecting tips.

1.4 Material: PC (Polycarbonate)

1.5 Color: White

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-37
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Fixed Premium Pipettes, PC (Non-Autoclavable)

1.6 Basic Features

1.6.1 Made up of highly Virgin thermoplastic polycarbonate

1.6.2 CE certified

1.6.3 Chemical resistance: see website (www.abdoslifesciences.com)

1.7. Pipettes

Code	Capacity	No. Per Case
E10001	5 µl	1
E10002	10 µl	1
E10003	20 µl	1
E10004	25 µl	1
E10005	50 µl	1
E10006	100 µl	1
E10007	200 µl	1
E10008	250 µl	1
E10009	500 µl	1
E10010	1000 µl	1

1.8 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-37
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Fixed Premium Pipettes, PC (Non-Autoclavable)

Warranty: 12 months from the date of Installation or 15 months from the date of Invoicing

Instruction for Use

Product can be used for food contact also.
Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimen.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure.
Consult local fire codes prior to storage of flammable liquids.



TECHNICAL DATA SHEET

Doc No – Ab/TDS/LH-37
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Fixed Premium Pipettes, PC (Non-Autoclavable)

Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

