

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-199
Rev No: 06 Date- 02/2022-23

Abdos Labtech Pvt. Ltd.
Safety Goggles, PC

Valid for Catalog No. U40101

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Safety Goggles, PC

1.1.1 Polycarbonate lens with side shield.

1.1.2 Impact resistant and lightweight frames for comfortable eyewear.

1.1.3 99.9% UVA / UVB protection.

1.2 Design: Safety Goggles

1.3 Performance Specifications:

Code	No. Per Case
U40101	2

1.4 Material: PC (Polycarbonate) confirm US FDA 21CFR

1.5 Color: Transparent Natural

1.6 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.7 Other Information

1.7.1 Excellent clarity and visibility.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-199
Rev No: 06 Date- 02/2022-23

Abdos Labtech Pvt. Ltd. Safety Goggles, PC

1.7.2 Moderate chemical stability.

1.7.3 Temperature range from -5°C to +55°C.

1.8 Chemical resistance: See website (www.abdoslifesciences.com)

1.9 Manufactured in Aseptic Environment

2.0 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure

Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-199
Rev No: 06 Date- 02/2022-23

Abdos Labtech Pvt. Ltd.
Safety Goggles, PC

Marketed By :

Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



Electronic IFU



Scan to visit website

Quality Manager

Abdos Labtech Private Limited

Employee Code 30110

