

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-222
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Specialty Indicator Tape for Steam Autoclave

Valid for Catalog No. U20802

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description: Specialty Indicator Tape for Steam Autoclave

1.1.1 These autoclavable tapes adhere to almost any dry surface-glass, plastic, metal, paper and rubber.

1.1.2 Versatile tape indicates completion of normal steam sterilization process.

1.1.3 All specialty steam indicator tape accepts, and written notations may be used to seal and identify almost all laboratory material before and after processing.

1.1.4 Change in the colour of the written matter from white to dark completion of autoclave cycle.

1.2 Design: Specialty Indicator Tape for Steam Autoclave

1.3 Performance Specifications:

Code	Size (Inch)	No. Per Case
U20802	1 x 500	1 Roll

1.4 Color: White

1.5 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-222
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Specialty Indicator Tape for Steam Autoclave

1.6 Other information:

1.6.1 **Chemical resistance:** See website (www.abdoslifesciences.com)

1.6.2 Manufactured in Aseptic Environment

1.6.3 **Intended use:** Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Instruction for Use

Storage temp. Max 30°C
Before using the product check crackness etc.
Check all the information printed on pouches or boxes.
Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.
DO NOT place any plastic labware in a flame.
DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure.
Consult local fire codes prior to storage of flammable liquids.

Marketed By :

Abdos Labtech Private Limited

Plot No. 1 Lakshwari Bhagwanpur Roorkee 247667 (UK) India
For Complaints/Feedback Contact Customer Care at above address or
Ph.: +91-011-41081897 or Email us at labtech@abdosindia.com



TECHNICAL DATA SHEET

Doc No – Ab/TDS/UTL-222
Rev No: 07 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Specialty Indicator Tape for Steam Autoclave

Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

