

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ MRB-129

Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Reversible Rack with Cover, PP

Valid for Catalog No. P20206, P20207, P20206A, P20207A

SRN (Single Registration Number): IN-MF-000022272

1. Description/Specification

1.1 Description:

Reusable Polypropylene (PP) Micro centrifuge Tube Rack places for MCT 1.5 ml, MCT 2.0 ml & MCT 0.5 ml

1.2 Design:

Two Step Rack (Reverse) for easy placement and working with micro centrifuge Tubes.

Both side contain 1.5 ml/2.0 ml & 0.5 ml place for easy transfer.

Compatible with 1.5 ml, 1.5 ml and 2.0 ml Centrifuge tubes

Each hole has molded in numbers for easy sample identification.

Holds 1.5/2.0ml micro centrifuge tubes on one side with smaller positions for the 0.5ml tubes on another side.

16 x 5 Array on the 80 Place rack for 0.5, 1.5- & 2.0-ml micro centrifuge tubes and
12 x 8 array on the 96 place rack for 0.5, 1.5 ml & 2.0 ml micro centrifuge tubes.

Marked Alphabetical in Y- axis and Numerical in X- axis for easy discrimination.

1.3 Performance Specifications:

1.3.1 Capacity:

P20206 – 96 Place (12x8 Array) for 1.5/2.0. ml MCT & 6 Place for MCT 0.5 ml Reverse side.

96 Place for 0.5 ml MCT & 6 Place for MCT 1.5/2.0 ml Reverse side.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ MRB-129

Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd. Reversible Rack with Cover, PP

P20207 – 80 Place (16x5 Array) for 1.5/2.0 ml MCT & 7 Place for MCT 0.5 ml Reverse side.

80 Place for 0.5 ml MCT & 8 Place for MCT 1.5/2.0 ml Reverse side.

1.3.2 Universal fit for different makes of MCT.

1.4 Material:

1.4.1 Polypropylene (PP), Hydrophobic, confirm US FDA 21CFR.

1.4.2 The product complies with USP class VI. Free from Heavy metals, it is safe to use in contact with food stuff, Pharmaceuticals & drinking water.

1.4.3 Product is manufactured without the use of Raw Materials of animal origin. Therefore BSE/TSE issue does not concern this product to the best of our knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color: Green, Yellow & Blue

1.6 Sterilization: NA

1.7 Quality control:

Products are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.8 Other Information: Reusable can be autoclaved up to 10 autoclave cycle.

1.9. Basic Features:

1.9.1 Made up of High pure virgin Polypropylene (PP)

1.9.2 Temperature range: -80⁰C to +121⁰C

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ MRB-129

Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Reversible Rack with Cover, PP

1.9.3 Autoclavability: Yes, Autoclavable at 15 psi and 121⁰C for 20 minutes

1.9.4 Chemical resistance: See website (www.abdoslifesciences.com)

1.9.5 Manufactured in Aseptic Environment

1.9.6 Dimension: 220 x 150 x 50 mm (**P20206**) & 240 x 85 x 50 mm (**P20207**)

1.9.7 Intended use: Widely used in scientific studies of molecular biology, clinical chemistry, biochemistry and any type of scientific laboratory.

Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121⁰C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.

Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimens.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure


Consult local fire codes prior to storage of flammable liquids.

TECHNICAL DATA SHEET

Doc No – Ab/TDS/ MRB-129

Rev No: 08 Date- 05/2024-25

Abdos Labtech Pvt. Ltd.
Reversible Rack with Cover, PP



Quality Manager
Abdos Labtech Private Limited
Employee Code 30110

